MICRO 360™

MULTI-PURPOSE MICROSCOPE

INSTRUCTION MANUAL

MODEL #44125

INTRODUCTION

Congratulations on your Celestron microscope purchase. Your new Micro360 Dual Purpose Microscope is a precision optical instrument made of high quality materials to ensure durability and long life. It is designed to give you a lifetime of enjoyment with minimal maintenance.

Before attempting to use your Micro360, please read through these instructions thoroughly to familiarize yourself with its functions and operations to maximize your enjoyment and usage. Refer to the microscope diagram to locate parts discussed in this manual.

Micro360 provides powers from 40x to 640x. It is ideally suited for examining specimen slides of yeasts and molds, cultures, plant and animal parts, fibers, bacteria, etc. You can also examine small and thin objects

at low powers such as coins, rocks, insects, PC Boards and many other objects. Eight prepared slides are included to get you started. You can also explore the exciting microscopic world by creating your own specimen slides using the included blank slides and tools. An excellent book to read that describes various activities, projects and techniques you can enjoy with your new microscope is "The World of the Microscope" by Usborne publishing.

The final section provides simple care and maintenance tips for you to follow to ensure your microscope provides you with years of quality performance, usage and enjoyment.

Note: *This product was designed and intended for ages 13 and older.*

SPECIFICATIONS

| | |
|--------------------------------------|---|
| STAGE | Plain Stage with Clips – 90 mm x 97 mm (3.5"x 3.8") |
| WF 10X AND WF 16X EYEPIECES | Glass Optics |
| FOCUSER | Coarse Focus |
| OBJECTIVES | All Glass Optics – see magnification chart for powers |
| ILLUMINATOR – TOP AND BOTTOM LED | Uses 4 AA Alkaline batteries (user supplied) |
| HEAD | Monocular with 45° Incline and 360° Rotatable |
| NOSEPIECE | Triple with click stop |
| CONDENSER | N.A. 0.65 |
| WEIGHT / DIMENSIONS (WITH BATTERIES) | 3.3 lbs (1.5 kg), 5.25"x8.5"x11.5" (133x216x292 mm) |

STANDARD ACCESSORIES WITH YOUR MICROSCOPE

- WF 10x Eyepiece
- 4x, 10x, 40x (spring loaded) Objective Lenses
- Top Illuminator – LED with Adjustable Brightness
- Bottom Illuminator – LED with Adjustable Brightness
- Needle Probe and Eyedropper
- WF 16x Eyepiece
- Filter Wheel/Diaphragm – 6 Positions
- 8 Prepared Slides and 4 Blank Slides
- 4 Cover Slips
- Tweezers

MAGNIFICATION TABLE

Use the following table to determine magnification using your microscope's different eyepiece/objective lens combinations.

| OBJECTIVE LENS | 4X | 10X | 40X |
|-----------------|-----|------|------|
| WF 10x Eyepiece | 40x | 100x | 400x |
| WF 16x Eyepiece | 64x | 160x | 640x |

SETTING UP MICRO360

1. Remove Styrofoam container from carton.
2. Carefully remove microscope and accessories from container and set on a table, desk or other flat surface.
3. Remove bag covering the microscope.
4. Remove cap from eyepiece tube.
5. Install four AA batteries (user supplied) under the microscope's base (see images below).

6. Insert the WF 10x eyepiece in the eyepiece tube.

MICRO360 OPERATION

Before viewing specimens, please read these sections thoroughly regarding focusing, changing power (magnification), using the stage and adjusting illumination.

ROTATING THE HEAD

Your microscope's head can be rotated 360°. You can view from any position by moving the head (while holding the arm) to a desired location.

VIEWING A SPECIMEN

Carefully place a specimen slide under the stage clips and center the specimen. Images will be observed upside down and reversed right to left, with practice centering objects will become easier to accomplish. Prepared slides, blank slides and tools are included with your microscope to get you started (see images below).

You are now ready to focus and view a specimen. Be careful not to damage a specimen slide or valuable object. When using higher powers while focusing, make sure the objective lens does not hit the slide or object being viewed.

FOCUSING & CHANGING POWER (MAGNIFICATION)

Place specimen slide (or object) directly under the objective lens and use the focus knob to focus on the specimen. For smaller objects, place on a flat material for easier viewing.

1. Always start with lowest power (4x objective lens and WF 10x eyepiece). This objective lens and eyepiece combination will yield 40x – the lowest but most useful power for viewing large or small solid objects.
2. For higher powers, rotate the nosepiece to change the objective lens to 10x or 40x. Hold the microscope above the nosepiece with one hand and rotate the nosepiece with the other hand until it clicks in position. Be cautious not to let the objective lens touch the specimen slide or object when changing to higher powers. **Note: You should first turn the focus knob to lower the stage.**
3. You can replace the WF 10x eyepiece with the WF 16x eyepiece to obtain three additional magnifications, including the highest power possible (640x).

4. At higher powers, your views will be greatly magnified but somewhat darker. The most enjoyable views are seen at lower powers which have a wider field of view and brighter illumination.

ILLUMINATION

To achieve the sharpest and best views, proper illumination (lighting) must be used:

1. To power the illuminator(s) on, turn the brightness control switch (located on both sides of base).
2. The top illuminator is used for viewing at low power (4x –objective lens) as higher powers (10x and 40x objective lenses) will block some or all the light. If you need to use high power to observe solid objects, use a bright secondary light (desk lamp, etc.) for direct illumination.
3. The bottom illuminator is used mainly for viewing specimen slides where light shines up through a hole in the middle of the stage through the slide.
4. Having both illuminators on at the same time will provide enough light for thick and irregular specimens.

ADJUSTING THE LIGHTING

Specimens of different sizes, thickness and color variations will require different levels of illumination. You can adjust illumination by turning the brightness control switches as needed.

When viewing a specimen that is not transparent or dark in color, you may need to increase the amount of light to resolve certain features or details. Increase illumination by turning the brightness control switches all the way to their highest setting.

Optimum lighting is found by making adjustments and experimenting. Each specimen may require slightly different illumination as well as same specimens viewed under different magnifications.

When lighting is no longer needed, turn illuminators off to preserve battery life.

USING FILTERS AND DIAPHRAGM

Before using Micro360, check and make sure no filters are in the optical path. Viewing is normally conducted without colored filters. However, to bring out increased levels of detail, experiment by using different backlighting filter colors for very bright transparent specimens. To change lighting color, rotate the filter wheel to the desired color – Red (R), Green (G), Blue (B). Each color filter is centered once a faint click stop is heard or felt. You may need to adjust the focus knob slightly for best results.

DIAPHRAGM – Within the filter wheel are three diameter openings – (1) 1 mm, (3) 3 mm and (6) 6 mm, which limit the amount of light passing through to the specimen. These openings are part of the diaphragm which allows you to change diameter sizes to maximize contrast, brightness, etc. The default setting is “6”.


CARE, MAINTENANCE AND WARRANTY


Your Celestron microscope is a precision optical instrument and should be treated with care at all times. Follow these care and maintenance suggestions and your microscope will need very little maintenance throughout its lifetime.

- When you are finished using your microscope, remove any specimens left on the stage.
- Turn off top and bottom illuminators when you are finished using your microscope.
- Remove batteries if you will not use your microscope for an extended period of time.
- Store the microscope in a dry and clean place.
- Be very careful if using your microscope in direct sunlight to prevent damage to the microscope or your eyes.
- When moving your microscope, carry it by the “arm”.
- Clean outside surfaces with a moist cloth.
- Never clean optical surfaces with cloth or paper towels as they can scratch optical surfaces easily.
- Remove dust with a camel’s hair brush or use an air blower.
- To clean fingerprints off optical surfaces, use a lens cleaning agent and lens tissue available at most photo outlets. While cleaning, do not rub in circles as this may cause streaks and scratches to occur.
- Never disassemble or clean internal optical surfaces. Repair work should be performed by qualified factory technicians or other authorized repair facilities.
- Use care when handling glass specimen slides, as edges can be sharp.

WARRANTY

Your microscope has a two year limited warranty. Please see Celestron’s website for detailed information on all Celestron microscopes at www.celestron.com.

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All rights reserved. / Printed in China 01-12
Torrance, California 90503 U.S.A.